



Quality control standard metals for drinking water, high range, 500 mL

Artikel-Nr.: 2833649

CHF Preis: Kontakt
Lieferung innerhalb von 2 Wochen

Ready to use - no dilutions required. Mixed standard includes: Copper: 2.5 mg/L Cu, Iron: 1.5 mg/L Fe, Manganese: 5 mg/L Mn. 500 mL Bottle

Technische Daten

Beschreibung:	Metals Drinking Water (High Range) Mixed Parameter QC Standard
Konzentration:	Kupfer 2,5 mg/L Cu
Plattform:	Lösung