



Test kit color disc chlorine dioxide, model CLD-2

Artikel-Nr.: 2842800

CHF Preis: Kontakt

Lieferung innerhalb von 2 Wochen

Gefahrstoffe



Monitor ClO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Technische Daten

Anzahl Bestimmungen:	100
Kleinste Schrittweite:	0.1
Koffertyp:	D
Messbereich:	0 - 2.0 ppm ClO ₂
Methode:	Color Disc/DPD
Modell:	CLD-2
Parameter:	Chlorine Dioxide
Plattform:	Color Disc
Standard-Messbereich:	0 - 2.0 ppm ClO ₂
Test-Kit Modell:	CLD-2 ²