



Refill kit for 32313. CO2 calibration starter kit

Artikel-Nr.: 32314

VERALTETER ARTIKEL

Dieser Artikel ist nicht mehr verfügbar.

Reagents to prepare a reference solution to calibrate CO2 analyzers. High carbonation 3 to 8 g/kg. Refill kit for 32313.

Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to $\pm 0.5\%$. Each reagent kit contains a certificate of analysis.

Quick and Easy to Use

Technische Daten

Lagerbedingungen:	Bei Raumtemperatur an einem trockenen Ort lagern
Messbereich:	3 - 8 g/kg
Messgerät:	Orbisphere 3658 Tragbarer CO2-Analysator
Parameter:	CO2