



Sension+ 5064 Robust conductivity cell for field use, attractive priced

Artikel-Nr.:

LZW5064.97.0002

CHF Preis:

Kontakt

Kein Versanddatum angezeigt

Robust conductivity cell with Titanium electrodes and Polycarbonat body for Sension+ portable meters. Cell constant $k=1\text{ cm}^{-1}$. Built-in Pt1000 for automatic temperature compensation. 1 m cable, MP5 connector.

Measuring range $5\text{ }\mu\text{S/cm}$ - 50 mS/cm , $0 - 80\text{ }^{\circ}\text{C}$. Robust conductivity cell for general aqueous application with attractive price position.

Attractive price

Built-in Pt1000 for automatic temperature compensation.

Technische Daten

Abmessungen (D x L):	12 mm x 85 mm
Elektroden-Anschluss:	Sension+ MP5
Elektroden-Typ:	Portabel
Genauigkeit:	0.50% with EC5 Meter
Gewährleistung:	6 Monate
Kabellänge:	1 m
Länge:	85 mm
Material:	Sensor Body Outside: Polycarbonate; Inside: Glass, Sensor Cell: Titanium
Messbereich LF:	$5\text{ }\mu\text{S/cm}$ - 50 mS/cm
Parameter:	Conductivity
Sensortyp:	2-polig, Titan, $k = 1\text{ cm}^{-1}$
Temperaturbereich:	$0 - 80\text{ }^{\circ}\text{C}$
Thermistor:	Pt1000