



Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Artikel-Nr.: 2546140

CHF Preis: Kontakt

Kein Versanddatum angezeigt

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Technische Daten

Durchmesser:	90 mm
Einheiten:	100/pk
Material:	Glass Microfiber
Porengröße:	1.5 µm
Volumen / Packungsgröße:	100 St.