



Radiometer Analytical XR110 Referenzelektrode (Roter Stab, Masseverbindung, Schraubkappe)

Artikel-Nr.: B20B112

CHF Preis: Kontakt

Kein Versanddatum angezeigt

The XR110 Reference Electrode has a glass body, Red Rod Ag/AgCl reference technology system, refillable electrolyte (KCl 3M), ground joint, and porous pin junction. The XR110 reference probe is recommended for Voltammetry applications and Educational purposes. The XR110 is a screw cap electrode and requires an external connection cable for connection to a meter. (Replace Obsolete #B20B110)

Technische Daten

Abmessungen (D x L):	8 mm x 120 mm
Diaphragma:	Porous plug
Elektroden-Anschluss:	Schraubkappe S7
Elektroden-Typ:	Laboratory reference electrode
Elektrolytlösung:	KCl 3M
Gewährleistung:	12 Monate
Inhalt:	Radiometer Analytical XR110-Referenzelektrode, 30 mL KCl 3M-Nachfülllösung und Betriebsanleitung.
Material:	Glass
Temperaturbereich:	-10 to 100°C

Inhalt

Radiometer Analytical XR110-Referenzelektrode, 30 mL KCl 3M-Nachfülllösung und Betriebsanleitung.

Erforderliches Zubehör

- CL111 Elektrodenkabel FX/FIL/1m/B (Item A94L111)
- 1 m Kabel für S7-Schraubsteckkopf (Typ -9) mit Stecker: BNC (Item A94L114)